


Product Name: Recombinant Human LIPC Protein
Catalog #: PHH4788



Summary

Name	LIPC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11150
Host	HEK293
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PCGDMNSFSQ GLCLSCCKGR CNTLGYHVRQ EPRSKSKRLF LVTRAQSPFK VYHYQFKIQF INQTETPIQT TFTMSLLGTK EKMQKIPITL GKGIASNKTY SFLITLDVDI GELIMIKFKW ENSAVWANVW DTVQTIIPWS TGPRHSGLVL KTIRVKAGET QQRMTFCSEN TDDLLLRPTQ EKIFVKCEIK

image



Background

Alternative Names	Hepatic triacylglycerol lipase (HL) (Hepatic lipase) (EC 3.1.1.3) (Lipase member C) (Lysophospholipase) (EC 3.1.1.5) (Phospholipase A1) (EC 3.1.1.32)
--------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------

Note

Product Name: Recombinant Human LIPC Protein
Catalog #: PHH4788



For Research Use Only , Not for Diagnostic Use