


Product Name: Recombinant Human PTPRG Protein
Catalog #: PHH4806



Summary

Name	PTPRG_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P23470
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEREHEEDGE KDSEKKEKSG VTHAAEERNQ TEPSPTPSSP NRTAEGGHQT IPGHEQDHTA VPTDQTGGRR DAGPGLDPDM

image



Background

Alternative Names Receptor-type tyrosine-protein phosphatase gamma (Protein-tyrosine phosphatase gamma) (R-PTP-gamma) (EC 3.1.3.48)

Note

For Research Use Only , Not for Diagnostic Use