


Product Name: Recombinant Human SLC5A3 Protein
Catalog #: PHH4843



Summary

Name	SC5A3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P53794
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEVKRRYMLA SPDVTSILLT YNLSNTNSCN VSPKKEALKM LRNPTEDEDVP WPGFILGQTP ASVWYWCADQ VIVQRVLA AK NIAHAKGSTL

image



Background

Alternative Names	Sodium/myo-inositol cotransporter (Na ⁺)/myo-inositol cotransporter) (Sodium/myo-inositol transporter 1) (SMIT1) (Solute carrier family 5 member 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use