


Product Name: Recombinant Human SLC39A6 Protein
Catalog #: PHH4896



Summary

Name	S39A6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13433
Host	HEK293
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DHEHSDHER HSDHEHHSEH EHHSDHDHHS HHNHAASGKN KRKALCPDHD SDSSGKDPRN SQGKGahrPE HASGRRNVKD SVSASEVTST

image



Background

Alternative Names Zinc transporter ZIP6 (Estrogen-regulated protein LIV-1) (Solute carrier family 39 member 6) (Zrt- and Irt-like protein 6) (ZIP-6)

Note

For Research Use Only , Not for Diagnostic Use