



Summary

Name	TMC5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UXY8
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEDPDYPDYS GSNRTQGYL KTQGYPDVPG PLNNPDYPGT RSNPYSVASR TRPDYPGSLA EPNYPRSLSN PDYSGTRSNA YSAARTSPD HPTSLPEPDY SEFQSHPYHR ASSRQPDYPG SQRNPDFAGS SSSGNYAGSR THPDHFGSLE PDYPGAQSNS DHPGPRANLN HPGSRKNLEH TSFRINPYAD SLGKPDYPGA DIQPNSPFF GEPDYPSAED NQNLPTWRE PDYSDAENGH DYGSSETPKM TRGVLSRTSS IQPSFRHRSD DVGSLWGEN DYPEGIEMAS MEMANSYGH LPGAPGSGYV NPAYVGESGP VHAYGNPPLS ECDWHKSPQG

image



Background

Alternative Names Transmembrane channel-like protein 5

Product Name: Recombinant Human TMC5 Protein
Catalog #: PHH5022



Note

For Research Use Only , Not for Diagnostic Use