


Product Name: Recombinant Human DCHS2 Protein
Catalog #: PHH5023



Summary

Name	PCD23_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6V1P9
Host	HEK293
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ISIQVLDVND NPPNFSSLSY HTHVKESTPL GSHITVVSAN DRDTGSHAEI IYNIISGNEK GHFYLEENTG VLYLIKPLDY EKMTKFTLTV QASDAEKKHF SFAVVFVSVL

image



Background

Alternative Names Protocadherin-23 (Cadherin-27) (Cadherin-like protein CDHJ) (Cadherin-like protein VR8) (Protein dachsous homolog 2) (Protocadherin PCDHJ)

Note

For Research Use Only , Not for Diagnostic Use