


Summary

Name	MEGF8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z7M0
Host	HEK293
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSKGEPKKYS LDPEEIENWV TEGPSEDEAV CVNCQNNSYG EKCESCLQGY FLLDGKCTKC QCNGHADTCN EQDGTGCPCQ NNTETGTCQG SSPSDRRDCY KYQCAKRES FHGSPLGGQQ CYRLISVEQE CCLDPTSQTN

image



Background

Alternative Names	Multiple epidermal growth factor-like domains protein 8 (Multiple EGF-like domains protein 8) (Epidermal growth factor-like protein 4) (EGF-like protein 4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MEGF8 Protein
Catalog #: PHH5055

