


Product Name: Recombinant Human ADGRG6 Protein
Catalog #: PHH5057



Summary

Name	AGRG6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86SQ4
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CRVVLSPSG TFTSPCYPND YPNSQACMWT LRAPTGYIIQ ITFNDFDIEE APNCIYDSL LDNGESQTKF CGATAKGLSF NSSANEMHVS FSSDFSQKK GFNASYIRVA VSLRNQKVIL PQTSDAYQVS VAKSISPEL SAFTLCFEAT KVGHEDSDWT AFSYSNASFT QLLSFGKAKS GYFLSISDSK CLLNNALPVK EKEDIFAESF EQLCLVWNS LGSIGVNFKR NYETVPCDST ISKVIPGNGK LLLSNQNEI VSLKGDYINF RLWNFTMNAK ILSNLSCNVK GNVVDWQND WNIPNLALKA ESNLSCGSYL

image



Background

Alternative Names	Adhesion G-protein coupled receptor G6 (Developmentally regulated G-protein-coupled receptor) (G-protein coupled receptor 126) (Vascular
--------------------------	--

Product Name: Recombinant Human ADGRG6 Protein
Catalog #: PHH5057



inducible G protein-coupled receptor) [Cleaved into: ADGRG6 N-terminal fragment (ADGRG6-NTF); ADGRG6 C-terminal fragment (ADGRG6-CTF)]

Note

For Research Use Only , Not for Diagnostic Use