



Product Name: Recombinant Human SULF2 Protein
Catalog #: PHH5082



Summary

Name	SULF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8IWU5
Host	HEK293
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGLGDAAQPR NLTKRHWPGA PEDQDDKDGG DFGTGGLPD YSAANPIKVT HRCYILENDT VQCDLDLYKS LQAWKDHKLH IDHEIETLQN KIKNLREVRG HLKKKRPEEC DCHKISYHTQ HKGRLKHRGS SLHPFRKGLQ EKDKVWLLRE

image



Background

Alternative Names Extracellular sulfatase Sulf-2 (hSulf-2) (EC 3.1.6.-)

Note

For Research Use Only , Not for Diagnostic Use