


Product Name: Recombinant Human TMEM123 Protein
Catalog #: PHH5094



Summary

Name	PORIM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8N131
Host	HEK293
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGLPHNSSAN STETLQHVPS DHTNETSNST VKPPTSVASD SSNTTVTTMK PTAASNTTTP GMVSTNMTST

image



Background

Alternative Names	Porimin (Keratinocytes-associated transmembrane protein 3) (KCT-3) (Pro-oncosis receptor inducing membrane injury) (Transmembrane protein 123)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use