


Product Name: Recombinant Human GLG1 Protein
Catalog #: PHH5197



Summary

| | |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | GSLG1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-hFc & C-His |
| Accession # | Q92896 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 38 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | SELEGQVISC LKLRADQRL SSDCEDQIRI IIQESALDYR LDPQLQLHCS DEISLCAEE AAAQEQTGQV EECLKVNLK |

image



Background

Alternative Names Golgi apparatus protein 1 (CFR-1) (Cysteine-rich fibroblast growth factor receptor) (E-selectin ligand 1) (ESL-1) (Golgi sialoglycoprotein MG-160)

Note

For Research Use Only , Not for Diagnostic Use