


Product Name: Recombinant Human LRRC37B Protein
Catalog #: PHH5252



Summary

Name	LR37B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96QE4
Host	HEK293
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APSDSSGINL SGFGGDQLEI QLTEQLRSLI PNEDVRKFMS HVIRTLKMEC SETHVQGSCA KLMLRTGLLM KLLSEQQEAK ALNVEWDTDQ QKTNYINENM EQNEQKEQKS

image



Background

Alternative Names Leucine-rich repeat-containing protein 37B (C66 SLIT-like testicular protein)

Note

For Research Use Only , Not for Diagnostic Use