



Product Name: Recombinant Human MTNAP1 Protein
Catalog #: PHH5272



Summary

Name	CQ080_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9BSJ5
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VYQSKPATLP RAKKMKGPIK DLIKAKGKEL ETENEERNSK LVVDKPEQTV KTFPLPAVGL ERAATTKADK DIKNPIQPSF

image



Background

Alternative Names Uncharacterized protein C17orf80 (Cell migration-inducing gene 3 protein) (Human lung cancer oncogene 8 protein) (HLC-8)

Note

For Research Use Only , Not for Diagnostic Use