



Product Name: Recombinant Human ADGRE3 Protein
Catalog #: PHH5279



Summary

Name	AGRE3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9BY15
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ATTILRDVES KVLETALKDP EQKVLKIQND SVAIETQAIT DNCSEERKTF NLNVQMNSMD IRCSDIIQGD

image



Background

Alternative Names	Adhesion G protein-coupled receptor E3 (EGF-like module receptor 3) (EGF-like module-containing mucin-like hormone receptor-like 3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use