



Summary

Name	PC11Y_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9BZA8
Host	HEK293
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NYSYELVLP S TNP GTVVFQV IAVDNDTG MN AEV RYSIVGG NTRDLFAIDQ ETGNITLMEK CDVTDLGLHR VLVKANDLGQ PDSLFSVVIV NLFVNESVTN ATLINELVRK SIEAPVTPNT EIADVSSPTS

image



Background

Alternative Names	Protocadherin-11 Y-linked (Protocadherin-11) (Protocadherin on the Y chromosome) (PCDH-Y) (Protocadherin prostate cancer) (Protocadherin-PC) (Protocadherin-22)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PCDH11Y Protein
Catalog #: PHH5280

