


Product Name: Recombinant Human IMPG2 Protein
Catalog #: PHH5281



Summary

Name	IMPG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZV3
Host	HEK293
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IEEATKPAGE QIAEFSIHLL GKQYREELQD SSSFHHQHLE EEFISEVENA FTGLPGYKEI RVLEFRSPKE NDSGVDVYYA VTFNGEAINN TTWDLISLHS NKVENHGLVE LDDKPTVVYT ISNFRDYIAE TLQQNFLLGN

image



SDS



WB

Background

Alternative Names	Interphotoreceptor matrix proteoglycan 2 (Interphotoreceptor matrix proteoglycan of 200 kDa) (IPM 200) (Sialoprotein associated with cones and rods proteoglycan) (Spacrcan)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IMPG2 Protein
Catalog #: PHH5281

