


Product Name: Recombinant Human SEMA6A Protein
Catalog #: PHH5291



Summary

Name	SEM6A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2E6
Host	HEK293
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTFEQDIERG NTDGLGDCHN SFVALNGHSS SLLPSTTTSD STAQEGYESR GGMLDWKHLL DSPDSTDPLG AVSSHNHQDK KGVIRESYLK

image



Background

Alternative Names Semaphorin-6A (Semaphorin VIA) (Sema VIA) (Semaphorin-6A-1) (SEMA6A-1)

Note

For Research Use Only , Not for Diagnostic Use