



Product Name: Recombinant Human DNAJC17 Protein
Catalog #: PHH5354



Summary

Name	DJC17_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST & C-His
Accession #	Q9NVM6
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QAAERTQKLD EKRRKVKLDL EARERQAQAQ ESEEEEEERS TRTLEQEIER LREEGSRQLE EQQRLIREQI RQERDQRLRG KAENTEGQGT PKLKLKWKCK

image



Background

Alternative Names DnaJ homolog subfamily C member 17

Note

For Research Use Only , Not for Diagnostic Use