



Product Name: Recombinant Human SLC12A6 Protein
Catalog #: PHH5398



Summary

Name	S12A6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHW9
Host	HEK293
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TLSSRHIDVC SKTKEINNMT VPSKLWGFCC NSSQFFNATC DEYFVHNNVT SIQGIPGLAS GIITENLWSN YLPKGEIIEK PSAKSSDVLG

image



Background

Alternative Names	Solute carrier family 12 member 6 (Electroneutral potassium-chloride cotransporter 3) (K-Cl cotransporter 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use