


Product Name: Recombinant Human ZAN Protein
Catalog #: PHH5449




Summary

Name	ZAN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9Y493
Host	HEK293
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTTPTEETTT SMEEPVIPTE KPSIPTEKPS IPTEKPTISM EETIISTEKP TISPEKPTIP TEKPTIPTEK STISPEKPTT PTEKPTIPTE KPTISPEKPT TPTEKPTISP EKLTIPTEKP TIPTEKPTIP TEKPTISTEE PTPTEETTI STEKPSIPME

image



SDS



WB

Background

Alternative Names Zonadhesin

Note

For Research Use Only , Not for Diagnostic Use