



Product Name: Recombinant Human PCDHA5 Protein
Catalog #: PHH5462



Summary

| | |
|------------------------------|--|
| Name | PCDA5_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-hFc & C-His |
| Accession # | Q9Y5H7 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 51 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | FPLEGASDLD IGANAQLRYR LNPNEYFDLD VKTNEEETNF LELVLRKSLD REETQEHRLI VIATDGGKPE LTGTVQLLIN VLDANDNAPE FDKSIYNVRL LENAPSGTLV IKLNASDADE GINKEIVYFF SNLVLDDVKS KFIINSNTGE IKVNGELDYE DYNSEINID AMDKSTFPLS GHCKVVVKLL |

image



Background

Alternative Names Protocadherin alpha-5 (PCDH-alpha-5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PCDHA5 Protein
Catalog #: PHH5462

