


Summary

Name	S4A7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y6M7
Host	HEK293
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRSSQESQRQ APELLVSPAS DDIPTVVIHP PEEDLEAALK GEEQKNEENV DLTPGILASP QSAPGNLDNS KSGEIKNGS GGSRENSTVD FSKVDMNFMR

image



Background

Alternative Names	Sodium bicarbonate cotransporter 3 (Electroneutral Na/HCO ₃ cotransporter) (Sodium bicarbonate cotransporter 2) (Sodium bicarbonate cotransporter 2b) (Bicarbonate transporter) (Solute carrier family 4 member 7)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC4A7 Protein
Catalog #: PHH5470

