


## Summary

<b>Name</b>	S4A4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-hFc & C-His
<b>Accession #</b>	Q9Y6R1
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	YPINSNFKVG YNTLFSCTCV PDPANISIS NDTTLAPEYL PTMSSTDMYH NTTFDWAFLS KKECSKYGGN LVGNNCNFVP

## image



## Background

<b>Alternative Names</b>	Electrogenic sodium bicarbonate cotransporter 1 (Sodium bicarbonate cotransporter) (Na <sup>+</sup> )/HCO <sub>3</sub> <sup>-</sup> cotransporter) (Solute carrier family 4 member 4) (kNBC1)
--------------------------	---

## Note

For Research Use Only , Not for Diagnostic Use