


Product Name: Recombinant Human LAG3 Protein
Catalog #: PHH5481



Summary

Name	LAG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P18627
Host	HEK293
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVWAQEGAPA QLPCSTIPL QDLSLLRRAG VTWQHQPDSG PPAAAPGHPL APGPHPAAPS SWGPRRRYT VLSVGPGLR SGRLPLQPRV QLDERGRQRG DFSLWLRPAR RADAGEYRAA VHLRDRALSC RLRLRLGQAS MTASPPGSLR ASDWVILNCS FSRPDRPASV HWFRNRGQGR VPVRESPHHH LAESFLFPQ VSPMDSGPWG CILTYRDGFN VSIMYNLTVL GLEPPTPLTV YAGAGSRVGL PCRLPAGVGT RSFLTAKWTP PGGGPDLLVT GDNGDFTLRL EDVSAQAQAGT YTCHIHLEEQ QLNATVTLAI ITVTPKSFSG PGS LGKLLCE VTPVSGQERF VWSSLDTPSQ RSFSGPWLEA QEAQLSQPW QCQLYQGERL LGAAVYFTEL SSPGAQRSGR APGALPAGHL

image



Background

Product Name: Recombinant Human LAG3 Protein
Catalog #: PHH5481



Alternative Names

Lymphocyte activation gene 3 protein (LAG-3) (CD antigen CD223) [Cleaved into: Secreted lymphocyte activation gene 3 protein (sLAG-3)]

Note

For Research Use Only , Not for Diagnostic Use