



Product Name: Recombinant Human GAS6 Protein
Catalog #: PHH5487



Summary

Name	GAS6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14393
Host	HEK293
Species	Human
Molecular Mass	73 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLRPRQRRAF QVFEEAKQGH LERECVEELC SREEAREVFE NDPETDYFYP RYLDCINKYG SPYTKNSGFA TCVQNLPDQC TPNPCDRKGT QACQDLMGNF FCLCKAGWGG RLCDKDVNEC SQENGGCLQI CHNKPFSFHC SCHSGFELSS DGRTCQDIDE CADSEACGEA RCKNLPGSYS CLCDEGFAYS SQEKACRDVD ECLQGRCEQV CVNSPGSYTC HCDGRGGLKL SQDMDTCEDI LPCVPFSVAK SVKSLYLGRM FSGTPVIRLR FKRLQPTRLV AEFDFRTFDP EGILLFAGGH QDSTWIVLAL RAGRLELQLR YNGVGRVTSS GPVINHGMWQ TISVEELARN LVIKVNDAV MKIAVAGDLF QPERGLYHLN LTVGGIPFHE KDLVQPINPR LDGCMRSWNW LNGEDTTIQE TVKVNTRMQC FSVTERGSFY PGSGFAFYSL DYMRTPLDVG TESTWEVEVV AHIRPAADTG VLFALWAPDL RAVPLSVALV DYHSTKKLKK QLVVLAVEHT ALALMEIKVC DGQEHVVTVS LRDGEATLEV DGTRGQSEVS AAQLQERLAV LERHLRSPVL TFAGGLPDVP VTSAPVTAFY RGCMTLEVNR RLLDLDEAAY KHS DITAHSC

image



Product Name: Recombinant Human GAS6 Protein
Catalog #: PHH5487



Background

Alternative Names

Growth arrest-specific protein 6 (GAS-6) (AXL receptor tyrosine kinase ligand)

Note

For Research Use Only , Not for Diagnostic Use