



Product Name: Recombinant Human SELL Protein
Catalog #: PHH5489



Summary

Name	LYAM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P14151
Host	HEK293
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAHHGTDCTW YHYSEKPMNW QRARRFCRDN YTDLVAIQNK AEIEYLEKTL PFSRSYYWIG IRKIGGIWTW VGTNKSLEEE AENWGDGEPN NKKNKEDCWE IYIKRNKDAG KWNDDACHKL KAALCYTASC QPWSCSGHGE CUEIINNYTC NCDVGYGYPQ CQFVIQCEPL EAPELGTMDG THPLGNFSFS SQCAFSCSEG TNLTGIEETT CGPFGNWSSP EPTCQVIQCE PLSAPDLGIM NCSHPLASFS FTSACTFICS EGTELIGKKK TICSSGIWS NPSPICQKLD KSFSMIKEGD

image



Background

Alternative Names	L-selectin (CD62 antigen-like family member L) (Leukocyte adhesion molecule 1) (LAM-1) (Leukocyte surface antigen Leu-8) (Leukocyte-endothelial cell adhesion molecule 1) (LECAM1) (Lymph node homing
--------------------------	---

Product Name: Recombinant Human SELL Protein
Catalog #: PHH5489



receptor) (TQ1) (gp90-MEL) (CD antigen CD62L)

Note

For Research Use Only , Not for Diagnostic Use