



## Summary

<b>Name</b>	GFRA3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O60609
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	DPLPTESRLM NSCLQARRKC QADPTCSAAY HHLDSCTSSI STPLPSEEPS VPADCLEAAQ QLRNSSLIGC MCHRRMKNQV ACLDIYWTVH RARSLGNYEL DVSPYEDTVT SKPWKMNLSK LNMLKPDSDL CLKFAMLCTL NDKCDRLRKA YGEACSGPHC QRHVCLRQLL TFFEKAAEPH AQGLLLCPCA PNDRGCGERR RNTIAPNCAL PPVAPNCLEL RRLCFSDPLC RSRLVDFQTH CHPMDILGTC ATEQSRCLRA YLGLIGTAMT PNFVSNVNTS VALSCTCRGS GNLQEECEML EGFFSHNPCL TEAIAAKMRF HSQLFSQDWP HPTFAVMAHQ NEN

## image



## Background

<b>Alternative Names</b>	GDNF family receptor alpha-3 (GDNF receptor alpha-3) (GDNFR-alpha-3) (GFR-alpha-3)
--------------------------	--

**Product Name: Recombinant Human GFRA3 Protein**  
**Catalog #: PHH5497**

---



**Note**

For Research Use Only , Not for Diagnostic Use