


Product Name: Recombinant Human GFRA3 Protein
Catalog #: PHH5498



Summary

Name	GFRA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O60609
Host	HEK293
Species	Human
Molecular Mass	80 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DPLPTESRLM NSCLQARRKC QADPTCSAAY HHLDSCTSSI STPLPSEEPS VPADCLEAAQ QLRNSSLIGC MCHRRMKNQV ACLDIYWTVH RARSLGNYEL DVSPYEDTVT SKPWKMNLSK LNMLKPDSDL CLKFAMLCTL NDKCDRLRKA YGEACSGPHC QRHVCLRQLL TFFEKAAEPH AQGLLLCPCA PNDRGCGERR RNTIAPNCAL PPVAPNCLEL RRLCFSDPLC RSRLVDFQTH CHPMDILGTC ATEQSRCLRA YLGLIGTAMT PNFVSNVNTS VALSCTCRGS GNLQEECEML EGFFSHNPCL TEAIAAKMRF HSQLFSQDWP HPTFAVMAHQ NEN

image



Background

Alternative Names GDNF family receptor alpha-3 (GDNF receptor alpha-3) (GDNFR-alpha-3) (GFR-alpha-3)

Product Name: Recombinant Human GFRA3 Protein
Catalog #: PHH5498



Note

For Research Use Only , Not for Diagnostic Use