


Product Name: Recombinant Human WFIKKN1 Protein
Catalog #: PHH5505




Summary

Name	WFKN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q96NZ8
Host	HEK293
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGSDCDIWDG QPVCRCRDRC EKEPSFTCAS DGLTYYNRCY MDAEACLRGL HLHIVPCKHV LSWPPSSPGP PETTARPTPG AAPVPPALYS SPSPQAVQVG GTASLHCDVS GRPPPAVTWE KQSHQRENLI

image



SDS



WB

Background

Alternative Names	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP, follistatin, immunoglobulin, Kunitz and NTR domain-containing protein)
--------------------------	---

Note

Product Name: Recombinant Human WFIKKN1 Protein
Catalog #: PHH5505



For Research Use Only , Not for Diagnostic Use