



Product Name: Recombinant Human PLXNC1 Protein
Catalog #: PHH5516



Summary

Name	PLXC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O60486
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VESEVDTELE VKIQKENDNF NISKKDIEIT LFHGENGQLN CSFENITRNQ DLTILCKIK GIKTASTIAN SSKKVRVKLG NLELYVEQES

image



Background

Alternative Names Plexin-C1 (Virus-encoded semaphorin protein receptor) (CD antigen CD232)

Note

For Research Use Only , Not for Diagnostic Use