


Product Name: Recombinant Human RXFP3 Protein
Catalog #: PHH5539



Summary

Name	RL3R1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9NSD7
Host	HEK293
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TMNKAAGGDK LAELFSLVPD LLEAANTSGN ASLQLPDLWW ELGLELPDGA PPGHPPGSGG AESADTEARV

image



Background

Alternative Names	Relaxin-3 receptor 1 (RLN3 receptor 1) (G protein-coupled receptor SALPR) (G-protein coupled receptor GPCR135) (Relaxin family peptide receptor 3) (Somatostatin- and angiotensin-like peptide receptor)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use