



Product Name: Recombinant Human TFA5 Protein
Catalog #: PHH5547



Summary

Name	TFA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-SUMO
Accession #	Q7Z5A7
Host	HEK293
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLLIAYCSQL AAGTCEIVTL DRDSSQPRRT IARQTARCAC RKGQIAGTTR ARPACVDARI IKTKQWCDML

image



Background

Alternative Names Chemokine-like protein TFA5-5

Note

For Research Use Only , Not for Diagnostic Use