



Product Name: Recombinant Human TMEM130 Protein
Catalog #: PHH5557



Summary

Name	TM130_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N3G9
Host	HEK293
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YELNLTDDSP ATTGAVVTIS ASLVAKDNGS LALPADAHLY RFHWIHTPLV LTGKMEKGLS STIRVVGHPV GEFVSVVWVT AADCWMCQPV ARGFVVLPI EFLVGDLVVT QNTSLPWPSS YLTKTVLKVS FLLHDPSNFL KTALFLYSWD FGDGTQMVTE DSVVYYNYSI IGTFTVKLVK VAEWEEVEPD ATRAVKQKTG DFSASLKLQE TLRGIQVLGP TLIQTFQKMT VTLNFLGSPP LTVCWRLKPE CLPLEEGECH PVSVAAYN LTHTFRDPGD YCFSIRAENI ISKTHQYHKI

image



Background

Alternative Names	Transmembrane protein 130
--------------------------	---------------------------

Product Name: Recombinant Human TMEM130 Protein
Catalog #: PHH5557



Note

For Research Use Only , Not for Diagnostic Use