



Product Name: Recombinant Human TRA Protein
Catalog #: PHH5565



Summary

Name	TRAR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P0DTU3
Host	HEK293
Species	Human
Molecular Mass	68 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSVQEAETVT LSCTYDTSES DYILFWYKQP PSRQMILVIR QEAYKQQNAT ENRFSVNFQK AAKSFLKIS DSQLGDAAMY FCAYRSVAVNA RLMFGDGTQL VVKPNIQNPD PAVYQLRDSK SSDKSVCLFT DFDSQTNVSQ SKDSVYITD KTVLDMRSMD FKSNSAVAWS NKSDFACANA FNNSIIPEDT FFPSPSSCD VKLVEKSFET DTNLNFQNL

image



Background

Alternative Names T cell receptor alpha chain MC.7.G5 (MC.7.G5 TRA) (TR alpha chain TRAV38-2DV8*01J31*01C*01)

Note

Product Name: Recombinant Human TRA Protein
Catalog #: PHH5565



For Research Use Only , Not for Diagnostic Use