



Product Name: Recombinant Human TMIGD1 Protein
Catalog #: PHH5570



Summary

Name	TMIG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UXZ0
Host	HEK293
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LTVNGKTENY ILDTTPGSQA SLICAVQNHT REEELLWYRE EGRVDLKSGN KINSSSVCVS SISENDNGIS FTCRLGRDQS VSVSVVLNVT FPPLLSGNDF QTVEEGSNVK LVCNVKANPQ AQMMWYKNSS LLDLEKSRHQ IQQTSESFQL SITKVEKPDN GTYSCIAKSS LKTESLDFHL IVKDKTVGVP

image



Background

Alternative Names Transmembrane and immunoglobulin domain-containing protein 1

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TMIGD1 Protein
Catalog #: PHH5570

