

Product Name: Recombinant Human TRIL Protein
Catalog #: PHH5571



Summary

Name	TRIL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L0X0
Host	HEK293
Species	Human
Molecular Mass	9 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQQPSPSVAA AAGPAPQSLD LHKKPQRGRP TRADPALAEP TPTASPGSAP SPAGDPWQRA TKHRLGTEHQ

image



Background

Alternative Names	TLR4 interactor with leucine rich repeats (Leucine-rich repeat-containing protein KIAA0644)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use