



Product Name: Recombinant Human TRABD2A Protein
Catalog #: PHH5574



Summary

Name	TIK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q86V40
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KGKSKKTSTR PTLSTIFAPK VPTLEVPAPE AVSSGHSTLP PLVSRPGSAD TPSEAEQRFR KKRRRSQRRP RLRQFSDLWV RLEESDIVPQ

image



Background

Alternative Names Metalloprotease TIK1 (EC 3.4.-.-) (TRAB domain-containing protein 2A)

Note

For Research Use Only , Not for Diagnostic Use