



Product Name: Recombinant Human PRADC1 Protein
Catalog #: PHH5578



Summary

Name	PADC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BSG0
Host	HEK293
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGFRIHDYLY FQVLSPGDIR YIFTATPAKD FGGIFHTRYE QIHLVPAEPP EACGELSNGF FIQDQIALVE RGGCSFLSKT RVVQEHGGRA VIISDNAVDN DSFYVEMIQD STQRTADIPA LFLGRDGYM IRRSLEQHGL PWAIISIPVN VTSIPTFELL QPPWTFW

image



Background

Alternative Names	Protease-associated domain-containing protein 1 (Protease-associated domain-containing protein of 21 kDa) (hPAP21)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PRADC1 Protein
Catalog #: PHH5578

