



Product Name: Recombinant Human PPY Protein
Catalog #: PHH5582



Summary

Name	PAHO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P01298
Host	HEK293
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLEPVYPGDN ATPEQMAQYA ADLRRYINML TRPRYGKRHK EDTLAFSEWG SPHAAVPREL

image



Background

Alternative Names	Pancreatic prohormone (Pancreatic polypeptide) (PP) (Obinepitide) [Cleaved into: Pancreatic hormone (PH); Pancreatic icosapeptide (PI)]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use