



Product Name: Recombinant Human TSPAN12 Protein
Catalog #: PHH5594



Summary

Name	TSN12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O95859
Host	HEK293
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPVQWSDMVT LKARMTNYGL PRYRWLTHAW NFFQREFKCC GVVYFTDWLE MTEMDWPPDS CCVREFPGCS KQAHQEDLSD LYQEGCGKKM YSFLRGTKQL

image



Background

Alternative Names	Tetraspanin-12 (Tspan-12) (Tetraspan NET-2) (Transmembrane 4 superfamily member 12)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use