



Product Name: Recombinant Human SLC6A19 Protein
Catalog #: PHH5605



Summary

Name	S6A19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q695T7
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRYDDCFSTN ILTLINGFDL PEGNVTQENF VDMQQRCNAS DPAAYAQLVF QTC DINAF LS EAVEGTGLAF IVFTEAITKM

image



Background

Alternative Names	Sodium-dependent neutral amino acid transporter B(0)AT1 (Solute carrier family 6 member 19) (System B(0) neutral amino acid transporter AT1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use