

Summary

Name	S22AP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q6T423
Host	HEK293
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAFILDHRCW VHILDNDTIP DNDPGTLSQD ALLRISIPFD SNLRPEKCRR FVHPQWKLIH LNGTFPNTSE PDTEPCVDGW VYDQSSFST IVTKWDLVCE

image



Background

Alternative Names Solute carrier family 22 member 25 (Organic anion transporter UST6)

Note

For Research Use Only , Not for Diagnostic Use