


Product Name: Recombinant Human CRH Protein
Catalog #: PHH5623



Summary

Name	CRF_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-SUMO
Accession #	P06850
Host	HEK293
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQPLDFFQPP PQSEQPQQPQ ARPVLLRMGE EYFLRLGNLN KSPAAPLSPA SSLLAGGSGS RPSPEQATAN FFRVLLQQLL LPRRSLDSPA ALAERGARNA LGGHQEAPER ERRSEEPPIS LDLTFHLLRE VLEMARAEQL AQAHSNRKL

image



Background

Alternative Names	Corticoliberin (Corticotropin-releasing factor) (CRF) (Corticotropin-releasing hormone)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use