


Product Name: Recombinant Human MDGA2 Protein
Catalog #: PHH5634




Summary

Name	MDGA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q7Z553
Host	HEK293
Species	Human
Molecular Mass	60 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YSYSLQWTQM NPDAVDRIVA YRLGIRQAGQ QRWWEQEIKI NGNIQKGELI TYNLTELIKP EAYEVRTPL TKFGECDSTI RVIKYSAPVN PHLREFHCGF EDGNICLFTQ DDTDNFDWTK QSTATRNTKY TPNTGPNADR SGSKEGFYMY

image



SDS



WB

Background

Alternative Names	MAM domain-containing glycosylphosphatidylinositol anchor protein 2 (MAM domain-containing protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use