



Product Name: Recombinant Human LRRC19 Protein
Catalog #: PHH5646



Summary

Name	LRC19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H756
Host	HEK293
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CNFTEKNYTL IPADIKKDVT ILDSYNQIT LNGTDTRVLQ TYFLLTELYL IENKVILHN NGFGNLSSLE ILNICRNSIY VIQQGAFLGL NKLKQLYLCQ NKIEQLNADV FVPLRSLKLL NLQGNLISYL DVPPLFHLEL ITLYGNLWNC SCSLFNLQNW LNTSNVTLEN ENITMCSYPN SLQSYNIKTV PHKAECHSKF PSSVTEDLYI HFQPISNSIF NSSNNLTRN SEHEPLGKSW

image



Background

Alternative Names Leucine-rich repeat-containing protein 19

Note

Product Name: Recombinant Human LRRC19 Protein
Catalog #: PHH5646



For Research Use Only , Not for Diagnostic Use