



Product Name: Recombinant Human LRP10 Protein
Catalog #: PHH5650



Summary

Name	LRP10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z4F1
Host	HEK293
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IIFPNHACED PPAVLLEVQG TLQRPLVRDS RTSPANCTWL ILGSKEQTVT IRFQKLHLAC GSERLTRSP LQPLISLCEA PPSPLQLPGG NVTITYSYAG ARAPMGQGFL LSYQDWLMC LQEEFQCLNH RCVSAVQRCD GVDACGDGSD EAGCSSDPFP GLTPRPVPSL PCNVTLEDFY GVFSPPGYTH LASVSHQSC HWLLDPHDGR RLAVRFTALD LGFGDAVHVY DGPGPPESSR LLRLTHFSN GKAVTVETLS GQAVVSYHTV AWSNGRGFNA TYHVRGYCLP WDRPCGLGSG LGAGEGLGER CYSEAQRCDG SWDCADGTDE EDCPGCPPGH FPCGAAGTSG ATACYLPADR CNYQTFCADG ADERRCRHCQ PGNFRRCRDEK CVYETWVCDG QPDCADGSDE WDCSYVLPRK

image



Background

Product Name: Recombinant Human LRP10 Protein
Catalog #: PHH5650



Alternative Names

Low-density lipoprotein receptor-related protein 10 (LRP-10)

Note

For Research Use Only , Not for Diagnostic Use