


Summary

Name	KERA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	O60938
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSVRQVYEVH DSDDWTIHDF ECPMECF CPP SFPTALYCEN RGLKEIPAIP SRIWYLYLQN NLIETIPEKP FENATQLRWI NLNKNKITNY GIEKGALSQ KLLFLFLED NELEEVPSPL PRSLEQLQLA RNKVSRI PQG TFSNLENLTL LDLQNNKLVD NAFQRDTFKG LKNLMQLNMA KNALRNMP PR LPANTMQLFL DNNSIEGIPE NYFNVIPKVA FLRLNHNKLS DEGLPSRGFD VSSILDQLS HNQLTKVPRI SAHLQHLHLD HNKIKSVNVS VICPSPSMLP AERDSFSYGP HLRYLRDGN EIKPPIMMAL MTCFRLQAV II

image



Background

Alternative Names Keratocan (KTN) (Keratan sulfate proteoglycan keratocan)

Product Name: Recombinant Human KERA Protein
Catalog #: PHH5665



Note

For Research Use Only , Not for Diagnostic Use