



Product Name: Recombinant Human KERA Protein
Catalog #: PHH5666



Summary

Name	KERA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O60938
Host	HEK293
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GIEKGALSQL KLLFLFLED NELEEVPSPL PRSLEQLQLA RNKVSRIPOG TFSNLENLTL LDQNNKLVD NAFQRDTFKG LKNLMQLNMA KNALRNMPPR

image



Background

Alternative Names Keratocan (KTN) (Keratan sulfate proteoglycan keratocan)

Note

For Research Use Only , Not for Diagnostic Use