


Product Name: Recombinant Human HLA-DQA2 Protein
Catalog #: PHH5675



Summary

Name	DQA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P01906
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VASYGVNFYQ SHGPSGQYTH EFDGDEEFYV DLETKETVWQ LPMFSKFISF DPQSALRNMA VGKHTLEFMM RQSNSTAATN EVPEVTVFSK FPVTLGQPNT LICLVDNIFP PVVNITWLSN GHSVTEGVSE TSFLSKSDHS FFKISYLTFL PSADEIYDCK VEHWGLDEPL LKHWEPEIPA

image



Background

Alternative Names	HLA class II histocompatibility antigen, DQ alpha 2 chain (DX alpha chain) (HLA class II histocompatibility antigen, DQ(6) alpha chain) (HLA-DQA1) (MHC class II DQA2)
--------------------------	---

Note

Product Name: Recombinant Human HLA-DQA2 Protein
Catalog #: PHH5675



For Research Use Only , Not for Diagnostic Use