



Product Name: Recombinant Human FBLN5 Protein
Catalog #: PHH5680



Summary

Name	FBLN5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9UBX5
Host	HEK293
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESNQCVDVDE CATDSHQCNP TQICINTEGG YTCSTDGYW LLEGQCLDID ECRYGYCQQL CANVPGSYSC TCNPGFTLNE DGRSCQDVNE CATENPCVQT CVNTYGSFIC RCDPGYELEE DGVHCSDMDE CSFSEFLCQH ECVNQPGTYF CSCPPGYILL

image



Background

Alternative Names Fibulin-5 (FIBL-5) (Developmental arteries and neural crest EGF-like protein) (Dance) (Urine p50 protein) (UP50)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FBLN5 Protein
Catalog #: PHH5680

